

Handheld Microscopes



Mid-Valley
STEM-CTE HUB



www.midvalleystem.org
midvalleystemctehub@linnbenton.edu
Linn-Benton Community College
Albany Campus - CC-212



Handheld Microscopes

The Handheld Microscope Kit provides students with the opportunity to explore the microscopic world through hands-on investigation. Featuring pocket microscopes, magnifying glasses, and prepared slides, this kit enhances observation skills, scientific inquiry, and curiosity about the natural world. Students can examine everyday objects, biological specimens, and their own discoveries up close, fostering an appreciation for detail and scientific exploration.



Grade Level

3rd - 12th grades

Group Size

groups of 2 students per microscope

Time Duration

30 minutes - 1 hour

Content of Kits

Components

- 15x pocket microscopes
- 15x magnifying glasses
- Prepared microscope slides
- Empty Slides
- Cover glass

Consumables

- AAA Batteries



Usage

Getting Started

1. **Introduce Magnification** – Discuss how magnifying glasses and microscopes work and explain their role in scientific discovery.
 2. **Practice Focusing** – Demonstrate how to adjust the focus on a pocket microscope for clear viewing of different objects.
 3. **Explore Prepared Slides** – Start with pre-made slides to help students practice viewing under controlled conditions.
 4. **Encourage Exploration** – Allow students to collect their own specimens and compare them under different magnification tools.
-

Storage

Keep microscopes in their protective cases and designated storage bin to prevent damage.

Troubleshooting

- **Battery Check** – Ensure that pocket microscopes have fresh batteries or are fully charged before use.
- **Lens Cleaning** – Use a soft microfiber cloth to clean lenses and avoid scratches that could blur the image.
- **Proper Handling** – Teach students to hold microscopes steadily and use both hands to avoid dropping or damaging them.



Activity Guide

Beginner

Introduction to Magnification

Students will start by using magnifying glasses to observe everyday objects such as leaves, fabric, or printed text. Have them describe what they see and compare it to the unaided eye, discussing how magnification reveals new details.

Intermediate

Exploring Prepared Slides

Using pocket microscopes, students will examine prepared microscope slides of biological specimens such as plant cells or insect parts. Students will record observations, sketch what they see, and discuss differences between various specimens.

Advanced

Investigating Everyday Materials

Students will collect their own samples, such as salt crystals, hair strands, or paper fibers, and prepare them for microscope observation. Have students document their findings, hypothesize how different materials are structured, and present their discoveries.

Extension Activities:

Microscopic Scavenger Hunt

Create a list of textures, patterns, and microscopic structures for students to find using their microscopes. Allow students to work individually or in teams to locate and document as many of the listed items as possible.

Science Journal Project

Students will keep a microscopy journal where they sketch and label their observations over a set period of time, building a collection of microscopic discoveries. Encourage them to research and write short descriptions about what they find. Adjust the subject of observation to accommodate for the amount of time students have with the equipment.



Learning Extensions

STEAM Connections: Science

Learning Objectives:

- Understand the principles of magnification and how microscopes improve scientific observation.
- Develop critical thinking and observation skills through hands-on examination of various specimens.
- Learn how to prepare and analyze samples under a microscope.
- Record and communicate scientific findings effectively.

Career Connections:

- **Biology & Life Sciences** – Introduces students to the study of cells, organisms, and ecosystems through microscopic observation.
- **Forensic Science** – Builds skills in analyzing microscopic evidence, such as fibers or fingerprints, for crime investigations.
- **Material Science** – Encourages exploration of the microscopic properties of different materials, leading to careers in engineering and product development.
- **Medical Research** – Provides a foundation for careers in pathology, microbiology, and biomedical research, where microscopes are essential tools.

Essential Employability Skills:

- Critical thinking
- Problem-solving
- Attention to detail
- Documentation
- Communication
- Observation





Resources and Accessibility

Safety Guidelines

- **Handle Glass Slides Carefully** – Teach students to hold slides by the edges to avoid fingerprints and breakage.
- **Avoid Direct Sunlight** – Remind students not to use magnifying glasses or microscopes to focus sunlight, as this can cause burns or fire hazards.
- **Keep Lenses Clean** – Wipe lenses with a soft cloth to maintain clear viewing and avoid damage.
- **Use in a Stable Area** – Ensure microscopes are placed on a stable surface to prevent accidental drops.

Accessibility

- **Tactile Learning** – Pair microscopic observations with hands-on models or textured materials for a multisensory experience.
- **Audio Descriptions** – Offer verbal descriptions of microscopic images for students with visual impairments.
- **Adjustable Stations** – Set up work areas at different heights to ensure accessibility for all students.

Library Catalog



Library Resources



Feedback

QR to feedback survey

